



Synchronize ease-of-use,
safety and interoperability

SmartFusion Series HSP-il infusion Pump



Easy to use

- Compact and light weight design allow easy moving and transferring.
- Caregivers can easily install or remove pumps thanks to rotatable pole clamp.



Modular design enables vertical stacking from 2-12 pumps without any tool.

- Automatic pump door



Enhanced safety

- With audible and visual alarm, caregivers can easily locate the alarm and warning.



Automatic anti-free-flow design



Up&down pressure sensors monitor pressure accurately



Double air bubble sensors detect air-in-line promptly

- Drug library: more than 1000 drugs in storage and recommended dose limits.
- History record: more than 30,000 records stored



simplified workflow

Wireless technology allows for connection to infusion management system HK-M1000 which features complete interoperability with HIS and CIS.

