

BLOOD AND INFUSION WARMER HSP2800

hospi

High performance dry heater for warming of:

- Infusion fluids
- Blood transfusion and blood products
- Blood return feeds / dialysis fluids
- Enteral/parenteral nutrition and rinsing solutions



Blood and Infusion Warmer HSP2800



FEATURES

Technology of control:

- Microprocessor integrated with temperature sensors provides high precise control

Safety system:

- Permanent running self-tests, 24 hours continuous operating
- Double independent over-heating protections and automatic cut off
- Visual and acoustic alarm for high temperature / low temperature/sensor fault

Dual channels:

- Each channel can work independently with own temp. setting and temp control
- Two channels can be combined to enlarge the warming capacity

Heating up to the patient:

- IV tube is completely wrapped in, no heat loss

User friendly interface:

- Big LED screen showing set temp., actual temp., heating time and fault situation
- Easy and quick to set up, ready to use within minutes.

Open system:

- Accepts standard IV tube, no special disposables needed

TECHNICAL SPECIFICATIONS

Model:	HSP2800
Temperature setting:	33 – 41 °C
Power supply:	ac 100-240V/50-60Hz
Power Consumption:	Max. 180VA
Type of protection against electric shock:	Class 1
Degree of protection against electric shock:	BF Applied part, Defibrillation-protected
Degree of protection against ingress of liquids:	IPX2
Temperature accuracy:	±1.0 °C
Overheat protection:	42 °C/43 °C
Low temperature alarm:	32 °C
Warming up time:	from 20 °C to 36 °C approx 2 min.
Operating mode:	Continuous
Dimension (W*D*H):	200x130x250mm
Net Weight:	2.3kg
Warming profile:	Length 1400mm, compatible to IV tube 3.5mm - 5.0mm OD. Contact Lateral Medical for other options

*Subject to technical modifications without notice

hospi

HOSPI SOLUTION LIMITED

WhatsApp: +86 139 6223 1292

Web: www.hospisolution.com

Email: info@hospisolution.com