


BB-100 Top Grade

- Microprocessor based servo controlled temperature system
- Control mode: Air mode
- Humidity is adjustable continuously
- Set temperature, air temperature, heating power are displayed separately by LED
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure
- >37°C temperature set function
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Single wall hood, 4 operating windows and 2 iris ports, 3 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet, cabinet base
- Optional function: X-ray cassette tray, LED phototherapy unit



BB-100 Standard



- Microprocessor based servo controlled temperature system
- Control mode: Air mode
- Humidity is adjustable in two levels
- Set temperature, air temperature, heating power are displayed separately by LED
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure
- >37°C temperature set function
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Single wall hood, 4 operating windows and 2 iris ports, 3 infusion seals
- Locking device for front door
- RS232 connector, 1 tray
- Optional function: X-ray cassette tray, LED phototherapy unit



Optional function for baby incubators

LED phototherapy unit

(Upside BL-30D, downside BL-100D)



Baby weight system

(BB-200 / BB-300 series)

Electric elevator

(BB-200 / BB-300 series)



Air&Oxygen Blender

(BB-200 / BB-300 series)

X-ray cassette tray

(Not available when use downside phototherapy)

