

BB-200 Luxurious

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Humidity control system
- Set temperature, air temperature, baby skin temperature, humidity, timer and heating power are displayed separately by LED
- Integral timer to record therapy time
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure, low water level
- >37°C temperature set function
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Single wall hood, 4 operating windows and 2 iris ports, 5 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet, 1 tray
- Optional function:
X-ray cassette tray; Baby weight system;
Electric elevator; LED phototherapy unit;
Double wall hood; Air and oxygen blender



Hospisolution Limited
Tel.: + 86 13912963835
Web: www.hospisolution.com
Email: info@hospisolution.com



BB-200 Top Grade

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Humidity display system
- Set temperature, air temperature, baby skin temperature, humidity, timer and heating power are displayed separately by LED
- Integral timer to record therapy time
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure
- >37°C temperature set function
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Single wall hood, 4 operating windows and 2 iris ports, 5 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet, 1 tray
- Optional function: X-ray cassette tray; Baby weight system ;
Electric elevator; LED phototherapy unit; Air and oxygen blender, Double wall hood.



BB-200 Standard

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Humidity is adjustable continuously
- Set temperature, air temperature, baby skin temperature, timer and heating power are displayed separately by LED
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure
- >37°C temperature set function
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Single wall hood, 4 operating windows and 2 iris ports, 5 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet, 1 tray
- Optional function: X-ray cassette tray; Baby weight system ;
Electric elevator; LED phototherapy unit; Air and oxygen blender, Double wall hood.

