

BB-300 Plus

- Air and skin temperature mode servo controlled by micro-computer
- 8" color LCD display touchscreen with skin temperature graph
- Humidity control system, oxygen concentration control system, baby weight system
- Set temperature, air temperature, baby skin temperature, humidity, oxygen concentration, SpO₂ and heating power are displayed separately on LCD
- Self - testing function, multiple audible and visual alarms
- >37°C temperature set function
- Double skin sensors
- Triple protection for over temperature with automatically cut off function
- The inclination of infant bed is adjustable
- Double wall hood, 4 operating windows and 2 iris ports, 5 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet
- Optional:
Electric elevator, CCD camera,
LED phototherapy unit, Air and oxygen blender



New humidifier



New access port



BB-300 Top Grade

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Humidity control and oxygen concentration control system
- Set temperature, air temperature, baby skin temperature, humidity, oxygen concentration, timer and heating power are displayed separately by LED
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure, low water level, O₂ concentration deviation
- >37°C temperature set function
- Embedded integrated sensor module
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Double wall hood, 4 operating windows and 2 iris ports, 5 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet, 1 tray
- Optional function: X-ray cassette tray; Baby weight system ;
Electric elevator; LED phototherapy unit; Air and oxygen blender



BB-300 Standard

- Microprocessor based servo controlled temperature system
- Control modes: Air mode and baby skin mode
- Humidity control system, oxygen concentration display system
- Set temperature, air temperature, baby skin temperature, humidity, oxygen concentration, timer and heating power are displayed separately by LED
- Integral timer to record therapy time
- Self - testing function, multiple audible and visual alarms
- Alarm functions: Power failure, over temperature, temperature deviation, fan failure, temperature sensor failure, low water level,
- >37°C temperature set function
- Embedded integrated sensor module
- Triple protection for over temperature with separate cut off device, more safety system
- The inclination of infant bed is adjustable
- Double wall hood, 4 operating windows and 2 iris ports, 5 infusion seals
- Independent locking device for front door
- RS232 connector, Oxygen inlet, 1 tray
- Optional function: X-ray cassette tray; Baby weight system ;
Electric elevator; LED phototherapy unit;
Air and oxygen blender

