

INFANT PHOTOTHERAPY UNIT SERIES



BL-360D

- LED blue bulbs as the radiant source, irradiating uniformly
- 360° phototherapy radiation more effective
- Small in size, light and high irradiation
- Three phototherapy modes: Upside phototherapy and Downside phototherapy can be used separately, also can be used together
- Therapy time of up and down side phototherapy units, air temperature, skin temperature are displayed separately
- Automatically record each therapy time and total therapy time
- Power off memory, alarm silence key, locking key and self-test function
- Four alarm functions: power failure, probe failure, over temperature, fan failure
- Count down timer
- Two sides observation windows
- Air spring structure make the hood open and close easily
- The aluminum alloy is applied on the base, with 2 drawers and 1 tray



BL-200D

- Double sides phototherapy for more effective radiation: Upside phototherapy and downside phototherapy can be used separately, also can be used together
- Irradiance intensity of upside phototherapy is adjustable in 3 grades: Low, medium, high
- LCD screen of upside phototherapy unit displays each therapy time and total therapy time with countup and countdown timer
- Upside phototherapy unit adopts big blue LED bulbs, which is, efficient and long life
- Downside phototherapy unit adopts blue LED bulbs source and independent air cooling fan
- Protect boards of infant bed could be folded down for clinical use
- Upside phototherapy unit can be adjusted 360° horizontally and 180° vertically
- The height of upside phototherapy is adjustable
- Stand column could be adjusted 360°
- Anti-rusty aluminum alloy base
- With drawers and lockable castors

