

BL-50D

- LED bulbs as phototherapy radiant source
- Irradiance intensity is adjustable in 3 grades: low, medium, high
- LCD screen display therapy time and integral time separately
- Count up timer and count down timer for exact and convenient treatment
- Equal light distribution, high intensity
- Quiet, no noise of fan
- The angle of head and height are adjustable
- Four castor with brake, aluminium alloy and steel support base
- Long life span LED bulbs
- Using together with infant incubators, infant warmers, baby cradle



BL-60D

- LED bulbs as phototherapy radiant source
- Irradiance intensity is adjustable in 3 grades: low, medium, high
- LCD screen display therapy time and integral time separately
- Count up timer and count down timer for exact and convenient treatment
- Equal light distribution, high intensity
- Quiet, no noise of fan
- The angle of head and height are adjustable
- Long life span LED bulbs
- Using together with infant warmers



BL-30D

- LED bulbs as phototherapy radiant source
- Irradiance intensity is adjustable in 3 grades: low, medium, high
- LCD screen display therapy time and integral time separately
- Count up timer and count down timer for exact and convenient treatment
- Equal light distribution, high intensity
- Quiet, no noise of fan
- Long life span LED bulbs
- Using together with baby incubator, installed on the hood



BL-100D

- LED bulbs as phototherapy radiant source
- Button control switch
- Equal light distribution, High intensity
- Quiet, two cooling fan, Long life span
- Timer for total using time and each treatment time
- Count down timer
- Irradiance intensity is automatic adjustment, to prevent the bed temperature higher
- Adopts Aluminum PCB to help element heat dissipation
- Using together with infant incubators, installed under the baby bed

