

# INFANT RADIANT WARMER

## BN-200 Standard

- Microprocessor based servo controlled temperature system
- Control modes: pre-warm,manual and automatic
- Audible and visual alarm functions: Power failure, temperature deviation, temperature sensor failure, skin temperature failure, over temperature
- 8" LCD colorful touch screen
- Carbon steel heater more effective
- Special design Inclined heater head is more convenient for X-ray image
- Integral timer to record therapy time
- Bigger mattress,
- Skin temperature trend for 2 hours
- Electric infant bed step-less adjustable tilt  $\pm 12^\circ$
- Four transparent protector with two infusion seals
- Convenient X-ray cassette tray under infant bed
- Electric elevator
- LED observation light, oxygen supply system and LED phototherapy
- With RS-232 connector, 2 trays and 2 drawers
- Optional Parts:  
low pressure aspirator ,internal / external infant resuscitator (T piece), Air and oxygen blender, Gel mattress



## BN-100 Luxurious

- Microprocessor based servo controlled temperature system
- Control modes: pre-warm, manual and automatic
- Set temperature and skin temperature display separately
- Special design Inclined heater head is more convenient for X-rays image
- Carbon steel heater more effective
- Audible and visual alarm functions: Power failure, temperature deviation, temperature sensor failure, skin temperature sensor, over temperature
- APGAR timer to record therapy time
- Convenient X-ray cassette tray under infant bed
- LED observation light adjustable, brightness adjustable
- Four transparent protector with two infusion seals
- LED phototherapy unit, oxygen supply system
- Bigger infant bed, Electric infant bed step-less adjustable tilt  $\pm 12^\circ$
- Built in low pressure aspirator
- With RS-232 connector, 2 trays and 2 drawers
- Optional Parts: internal / external infant resuscitator (T piece) ,  
Air and oxygen blender ,Gel mattress, Electric elevator,  
Air temperature & count down timer, Infant head fixing unit,  
Infant oxygen hood.

