

Half Sized, Same Wised

hospi

SE-1200 Pro


Electrocardiograph




EDAN

SE-1200 Pro Electrocardiograph

Designed for the convenience of caregivers, the SE-1200 Pro enables easy-going workflows. The device will automatically record ECG once the signals are well sampled, and present a timely reminder when the lead wires are misconnected, or the pacemaker is detected.



Intellectual
Sampling



Pacemaker
Detection



Self-adapt
AC Filter



Serious
Diseases Hint



Reverse
Electrode Hint

Seamless Connectivity



Hospitals
EMR/HIS/PACs




Clinics
GDT-based

DICOM/HL7

GDT



SE-1515 Data Management System





Serious Diseases Notification

A red notification will identify the severe illnesses after the diagnosis.

Pacemaker Analysis

80 kHz sampling rate.
 Auto pacemaker detection.
 Pick up pace signals lowest from $30\mu\text{s}$, $500\mu\text{V}$.
 Independent pacemaker detection channel.




Averaged Template

The device identifies the averaged template's P, Q, S, and T points with the relevant measurement parameters. Doctors can relocate these points manually and recalculate the corresponding parameters.

Gesture Amplification & Measurement

Users can amplify any segment of waveforms and measure it on the screen with gesture operation.



Reports Comparison

Users can compare two reports on the screen to follow a patient's clinical condition.



SE-1200 Pro

Electrocardiograph