

Vein Finder VT30



**NEW
Model**

VT30 Parameter

Light Type	Near-infrared light
Best Projection Distance	200±20mm
Dimension	240x60x60mm(LxWxH)
Battery Working Time	4 hours
Operating Mode	Handheld or suspension

VT30 Optional Stand



Visualization Modes

Brightness

Adjustable brightness allows clinicians customizing the image brightness base on the environment light.


Colors

Switchable colors can fit different needs from the nurses.

Inverse

Clinicians can highlight the veins or highlight the surrounding tissue.

Application Images:



Application Department:

Pediatrics 	Oncology 	Radiology 	Emergency 	Geriatric 	Laboratory 	Surgery 	Outpatient 	Plastic 	Vascular
----------------	--------------	---------------	---------------	---------------	----------------	-------------	----------------	-------------	--------------