


Endless Love 

To Guard The New Life Safe



HSP-D60/D70/D80 Patient Monitor

hospi

Product configuration ("✓" represents standard function, while "●" stands for optional accessory.)

	ECG	SpO ₂	Respiration	Non-invasive blood pressure	Temperature	End-tidal carbon dioxide	Print	Screen size
D60	✓	✓	✓	✓	✓	●	●	8 inches
D70	✓	✓	✓	✓	✓	●	●	10 inches
D80	✓	✓	✓	✓	✓	●	●	12 inches



HD screen for clear display



Touch (optional) & shuttle dual operation for work efficiency improvement



Extra wide viewing angle for observation from different angles



Intelligent light control of the screen to automatically adapt to the light environment



Fanless design to create a quiet environment for the department



Easy to disassemble the plug-in box for easy functional expansion



360° alarm lamp, easy to observe alarm



Optional mounting brackets, convenient to combine with other medical devices



Large font interface

Display the measured values of each parameter in an eye-catching bold type, easier for medical staff to identify.

Accurate and reliable physiological parameter monitoring

■ SpO₂ monitoring

Non-invasive, continuous measurement of pulsating oximetry was used to realize the simultaneous monitoring of blood oxygen parameters in patients.

■ CO₂ monitoring (optional)

Provide two measurement methods: mainstream and external bypass.

■ ECG monitoring

Professional ECG algorithm, suitable for monitoring the 3-lead (optional configuration) and 5-lead (standard configuration) ECG waveform and heart rate.

■ NIBP monitoring

Providing dual overvoltage protection for hardware and software to safe and accurate blood pressure monitoring for children and adults over 3 years of age.



hospi

HOSPISOLUTION LIMITED

Whatsapp: +86 139 6223 1292

Web: www.hospisolution.com

Email: info@hospisolution.com