


hospi

hospi

www.hospisolution.com

HSP-C2118

Mobile Digital FPD C-arm System




hospi

HOSPISOLUTION LIMITED

Whatsapp: +86 139 6223 1292
Web: www.hospisolution.com
Email: info@hospisolution.com

TECH
CONTRIBUTES
HEALTH



Less is more



All-in-One Design

The integrated design of C-arm and workstation allows for compact frame structure, small footprint, flexible and light movement. It can be operated by a single person and is suitable to use in crowded operating rooms, saving room space and facilitating doctors' clinical operations.



More flexibility, More possibilities

Indifferent equilibrium design

The C-arm frame can be hovered at any angle in the unlocked state to ensure safe and stable operation of the equipment.

27-inch ultra-wide screen display

The large-size 4K LCD display creates clear and delicate images; the display can be rotated at a large angle to meet the needs of clinical observation of images from multiple angles.

Touch screen

Equipped with high-quality and large-size touch operation screen, high sensitivity, flexible rotation, convenient for clinicians to perform touch operation from all angles of the machine.

One-key switch

The machine can be started quickly and in time when urgently needed, which can improve work efficiency extremely.

10-minute standby power

The machine can realize the transition between multiple operating rooms without power connection, which improves the working efficiency of the machine.

Two-way laser positioning

Equipped with laser localizers at both the tube and flat panel detectors to meet different positioning requirements.



Sharp in every image, Care for each clinician

Dynamic flat panel detector

Equipped with a 9 inches * 9 inches flat panel detector, output grayscale 16 bit, high DQE, low noise, dynamic range promotion, superior imaging performance, clear images, and meet various clinical examination needs.

High frequency inverter technology

The ray waveform is a rectangular wave with low scattered rays, which prevents the generation of soft rays from the source, ensures the ray dose control, and effectively cares for the health of doctors and patients.

With low dose mode

The machine is set with low-dose exposure parameters to strictly control unnecessary radiation dose.

Virtual collimators

The virtual collimators design can realize the adjustment of the size of the exposure range to reduce the radiation caused by repeated clinical exposure in a wireless state. And automatically adjust the size and location of the image interest area, using the automatic mode is more accurate.

Pluggable grid

The pluggable grid complies with national industry standards and is used in pediatrics and other dose-sensitive applications. Easy to use, while ensuring clear images, it can effectively protect against radiation hazards.

Intelligent frequency conversion technology

Automatically adjust the image frame rate according to the body part and the radiation dose to reduce the radiation dose while ensuring the image quality.

4,096

Traditional image intensifier

65,536

Gray-scale degrees

Dynamic plate detector





Intelligent image processing technology

Intelligent real-time image processing technology automatically analyzes images and performs enhancement processing, reducing tedious operations when processing images in clinical applications. It provides high-quality images quickly and completes clinical work efficiently.

1-second imaging

The design of images coming out within 1 second can save doctors' time in the clinical application process, thereby improving work efficiency.

DAP

DAP can be directly displayed during use to protect the health of doctors and patients..

Wireless exposure function (Optional)

The optional wireless foot pedal can achieve compartment exposure.

Clinical Image



Lumbar vertebrae



Lateral lumbar



Hip orthosis



Lateral cervical



Knee joint



Lateral ankle