

YA-PS09B

Patient Transfer Stretcher



YA-PS09B PATIENT TRANSFER STRETCHER

The YA-PS09B stainless steel patient transfer stretcher adopts full 304 stainless steel integrated frame, lightweight yet sturdy structure, suitable for ward, emergency room, clinic and first-aid transfer scenarios. Manual multi-angle adjustment design, matched with full set of medical auxiliary accessories, realizes safe, convenient and stable patient transportation, meeting daily basic medical transfer and examination demands.

Main features

- High-quality 304 stainless steel structure, 1.0mm wall thickness of stainless steel round tubes
- Four-section stainless steel bed panel with restraint strap holders on head and leg sections
- Back section adjustable from 0-75°, controlled by foldable hand screw rod with chrome-plated crank
- Leg section adjustable from 0-45°, controlled by independent foldable hand screw rod with chrome-plated crank
- 5-position aluminum guardrail with stainless steel upright, foldable toward tail end and detachable
- Height-adjustable stainless steel infusion pole with fixing hole at head end
- 125mm silent full-brake casters with straight plug-in rod fixation
- Bottom frame equipped with stainless steel tray and oxygen cylinder fixing rack
- 60mm thick one-piece PU medical mattress
- Screw-fixed detachable support legs, easy to assemble

Main details



Bed Platform

Four-section stainless steel bed panel, with restraint strap holders on the head and leg sections. Equipped with 60mm thick one-piece PU medical mattress for patient comfort and easy cleaning.



Guardrail

5-position aluminum guardrail with stainless steel uprights, which can be folded flat toward the tail end and is detachable, providing reliable anti-fall protection during transfer.



Pull-out Instrument Tray

Pull-out instrument tray stores medical tools, dressings and emergency supplies, easy to take and clean.



Oxygen Cylinder Bracket

Stainless steel oxygen cylinder bracket stably fixes oxygen tank to avoid shaking during moving.



Recessed Handle

Dual recessed handles on the side separately control backrest and leg plate lifting, simple operation, medical staff can adjust patient posture quickly without extra tools.



Universal Casters with Brake

125mm silent full-brake casters, fixed by straight plug-in rods, with flexible steering and stable parking, ensuring quiet and safe movement in hospital environment.



Infusion pole

Height-adjustable stainless steel infusion pole, with reserved fixing hole at the head section, can be installed or removed according to actual use requirements.

Backrest & Leg Section Adjustment

Backrest and leg section are controlled by independent hand cranks, stepless angle adjustment. Backrest can be lifted vertically for sitting posture support, leg plate can be raised to lift lower limbs, convenient for infusion, wound care and patient rest. Multiple tilt postures can be switched freely as shown in different states.



Technical parameter

Model	YA-PS09B
Overall size	1900*750*700mm
Material	High-quality 304 stainless steel, round tube wall thickness 1.0mm
Back section adjustment	0-75°
Leg section adjustment	0-45°
Guardrail	5-position aluminum guardrail, stainless steel upright, foldable & detachable
Infusion pole	Stainless steel, height adjustable
Castor	Dia 125mm, silent full-brake caster
Mattress	60mm thick one-piece PU medical mattress

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 30 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia. Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval. Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US