

YA-PS10A

Patient Transfer Stretcher



*Better care
Better comfort!*

YA-PS10A

PATIENT TRANSFER STRETCHER


DESCRIPTION

Stainless steel frame, heavy-duty structure, and compact design, all the features make YA-PS10A a practical patient transfer stretcher an affordable cost for durable operation.

- FEATURES
- High-quality stainless steel material
 - Has back section adjustment
 - Stainless steel long-side rail
 - Large castors easy to move
 - Washable and disinfectable leather mattress
 - Detachable structure reduces transportation costs

STANDARD SPECIFICATION

Overall size	L1850*W650*H750 mm
Mattress thickness	5cm
Back section adjustment	0 ~ 75°±5°
Castor	Dia 125mm
Safe working load	150 KG



Back section adjustment

The back section is adjusted by rack from 0~60°




Side rails

Stainless steel long side rails that rotate and drop down, ensuring patient safety.



Waterproof mattress

A washable and disinfected leather mattress, waterproof and compliant with hospital standards



Wheel

Four castors with diagonal braking ,dia 125 mm



Knock down design

A detachable structure ensures stability while significantly reducing transportation costs

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US