

YA-PS10B

Patient Transfer Stretcher



YA-PS10B PATIENT TRANSFER STRETCHER

The YA-PS10B adopts powder-coated frame structure and comes equipped with a removable top table, standard adjustment functions and smooth casters. Cost-effective fully functional, it fulfills core demands for patient transfer and bedside care, while the under-bed storage compartment further boosts space utilization flexibility.

Main features

- Powder coating structure, easy to clean and durable
- Back section can be adjusted by ratchet
- Side rail can be raised and dropped down
- Mattress thickness 5cm
- Large storage shelf
- Four castors with diagonal lock
- Knock down packing

Main details



Back section adjustment

Back section adjusted by ratchet and can be stopped at any angle.



Side rail

Rotatable height-adjustable armrests with bayonet locking for firm fixation.



Storage space

Big storage shelf at bottom of the stretcher.



Top table

The top of the transport stretcher can be removed, making it easier to move patients onto stretchers and reducing damage during transport.



Castors

4 units silent castors with diagonal lock, easy to move.



Bumpers

Bumpers are on the four corners, which can reduce the impact of accidentally bumping into the wall during the movement.

Technical parameter

Model	YA-PS10B
Overall Size(LWH)	1850×600×800 mm(72.8"x 23.6" x 31.5")
Mattress thickness	5cm (2.0")
Packing	Knock down design
Safe Working Load	180 kg
Castor	Dia 125mm (4.9")

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia. Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval. Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US